Communicator 5

Setting up Email:

1. Open the quick menu

(touch and hold anywhere on the screen or press this button on the left side of the device. You can also right click if you are using a mouse. Or select Ctrl + M on the keyboard)





- 2. Select Caregiver Settings
- 3. Select Applications tab
- 4. Select the email button
- Select new and then choose your email provider*.

(If you email provider is not listed, you will choose the configure manually option, if you need help filling in the information, contact your email provider to get the correct server information)

- Enter your name, email address and password (same password as you sign into your email)
- 7. To access email you will need to have access to an email page set from your home page
- 8. To add an email page, open the quick menu
- 9. Select All Page Sets button
- 10. Select Text Communication
- 11. Select Long Distance Communication
- 12. Select an email page that best meets your needs
- 13. Use the functions available on the email page to manage, send and receive email messages.

*If using Gmail, you will need to enable POP:

- Sign in to Gmail.
- Click the gear in the top right.



- Click Forwarding and POP/IMAP.
- Select Enable POP for all mail or Enable POP for mail that arrives from now on.
- Choose what to do with your messages after your POP client or device receives them.
- Click Save Changes.

Although the radio button next to your selected option may not stay selected, Gmail displays the status next to 'POP Status:.'

When you enable POP, your settings display 'POP Status: POP is enabled.'

If your POP settings are turned off, the status will indicate that POP is disabled.

Once you've enabled POP in your Gmail settings, you need to configure in Communicator 5.

For more information watch the video at

Watch dynavoxtech.force.com/devices/apex/Communicator

Please note you will need a mytobiidynavox.com account to access the resources



