



# Tobii Dynavox I-110

## Getting Started Guide




# What's in the Box?




- 1 Device
- 2 Power Cable
- 3 Shoulder Strap
- 4 Getting Started Guide
- 5 Safety & Compliance, License, and Configuration Documents


## Getting to Know Your Device



Front



Sides



Back

- A Power Button
- B Volume Buttons
- C Home Button

- D Audio & Switch Ports
- E Power Connector & USB Ports

- F Stand
- G Mounting Location
- H Strap Connectors

Follow the three steps in this guide and your new I-110 device will be set up and ready to use.

# Step 1: Initial Setup

## Start the Device


- 1 Connect the power cable to the device and plug in to an outlet.
- 2 Press the Power button.



## Setup


When you turn on your device, you will see one of two screens. If your device has not been pre-configured, you will be prompted to complete the Windows Setup process first.

### Windows Setup



While this process does not require a Microsoft account, we recommend that you connect to or create a user's account and not a caregiver's account.

### Start Up Wizard: Pick Your App



Choose which app to automatically launch when the device is turned on.

# Step 2: Communication Software Configuration



**TD Snap**  
Symbol-Supported  
Communicators

## TD Snap Setup Wizard

Follow the prompts to create your new user, or restore a saved user file.



**Communicator 5**  
Text-Supported  
Communicators

## Communicator 5 Setup Wizard

Follow the prompts to create your new user.

## Step 3: Mount & Position

The device may be positioned using a mounting system, laying it flat on a surface, or propped on the integrated kickstand. Start by situating the user comfortably, then find the device position that gives them clear screen visibility and easy access to their selection method of choice. Always position the device to suit the user, not the other way around. It is expected that the device may need to be repositioned throughout the day. Attach the shoulder strap for secure transport.



## Access Method Tips for TD Snap



### Touch Screen Access Methods

- Touch, Touch Enter, and Touch Exit allow the user to navigate and make selections using the touch screen on the device.
- If swiping or selecting buttons is difficult, consider changing the access method to Touch Enter or Touch Exit.
- Keyguards are available for all pre-set grid size configurations.



### Mouse Dwell

- Navigate using a head mouse, traditional mouse, or joystick plugged into the USB port on your device.
- Use the mouse to move around the screen, then make selections using dwell or a switch.
- Standard mouse options may be difficult for some people. You might consider a specialty mouse such as a trackball, wireless mouse, or joystick.



### Scanning

- Choose the TD Scanning page set in the Setup Wizard. These pages have been designed to improve efficiency.
- If you did not select the Scanning page set in the setup process, create a new page set:
  - 1 Select **Edit**.
  - 2 Select **Page Set** tab.
  - 3 Select the current page set.
  - 4 Select **Create Page Set**.
  - 5 Select **Core First Scanning Page Set** and select **Next**.
  - 6 Select **Create**.
  - 7 Select **Core First Scanning** to load the Page Set.
- Download the Scanning Implementation Guide for information and ideas for teaching scanning skills.



**For individuals needing visual or auditory feedback, turn on Highlighting or Selection Feedback Sound.**

# We Are Here to Help



## myTobiiDynavox

Use your free account to Sync, share Page Sets, store backups, and more! Visit:

**[www.myTobiiDynavox.com](https://www.myTobiiDynavox.com)**



## Learning Hub

Tobii Dynavox Learning Hub:

**[learn.tobiidynavox.com](https://learn.tobiidynavox.com)**



## Training Materials

User Manual:

**<https://qrco.de/I-110manual>**

Communicator 5 Training Cards:

**<https://qrco.de/C5Cards>**

TD Snap Core First Training Cards:

**<https://qrco.de/TDSnapCards>**

TD Snap Text Training Cards:

**<https://qrco.de/TextCards>**

TD Snap Aphasia Training Cards:

**<https://qrco.de/AphasiaCards>**

TD Snap Scanning Implementation Guide:

**<https://qrco.de/ScanGuide>**



## Additional Support

Community:

**<https://qrco.de/TDFB>**

North America Technical Support:

**1-800-344-1778 ext. 1**

**For more  
information:**



**[qrco.de/I-110Support](https://qrco.de/I-110Support)**

**For more  
information about  
access methods:**



**[qrco.de/bbA7US](https://qrco.de/bbA7US)**

