





TD Pilot with TD Snap

Getting Started Guide



What's in the Box



1 Pre-assembled device

(iPad Pro 12.9, TD Pilot Base, Protective Case, ConnectIT/Readapt mount plate, Connection Cable USB-C to USB-C, Pre-Installed Battery)

2 Connection cable
Lightning – USB C

3 Charger with cable

4 Screwdriver

5 Getting Started Guide

6 Safety and Compliance document

7 Adjustable Bracket

8 TD Snap Training Cards

9 Warranty documents

10 iPad Box
(containing iPad charger)


11 Carry bag

Let's Get Started

This Getting Started Guide is designed for two people: the person who will use the device to communicate and a helper who can touch the screen and move the device. The sections with the hand icon must be completed by the helper. The sections with the eye icon should be completed by the person who will access the device using eye gaze. When you have completed the steps in this guide, the individual will be able to use eye gaze to communicate in the TD Snap app.



Power On



- 1 Connect the charger cable to the charging port on the side of the TD Pilot base, then plug the power cord into a socket.
- 2 Press the power button on the back of the TD Pilot to power it on.
- 3 Press the power button on the iPad to power it on.
- 4 Follow the on-screen prompts to set up the iPad.

When you have completed your iPad setup prompts, you'll see the Home screen containing your app icons.

TD Pilot Power Management

The TD Pilot must be powered on when connected to external power in order for the iPad to charge too.

When the TD Pilot is powered on, the eye tracker is always active, even when the iPad screen is locked. Because of this, the batteries will drain when the TD Pilot is turned on and not connected to external power. It is best to leave the TD Pilot turned on and plugged into the charger whenever you are not using it.

Battery status information is found in *TD CoPilot > Settings > Battery*.

Additional information about power management is available in the TD Pilot User's Manual.

Set Up Eye Gaze Access

Configure the display for gaze accuracy

1 Tap **Settings**.



- 2 On the left side, tap the **Home Screen & Dock** category.
- 3 On the right side, select **Use Large App Icons** to make the icons larger.
- 4 On the left side, tap the **Display & Brightness** category.
- 5 On the right, tap **Dark**.
- 6 On the right, scroll down and tap **Text Size**.
- 7 Move the Text Size slider all the way to the right.



Note

This increases the text size in all compatible applications.

- 8 Tap **Display & Brightness** on the left side to return to the Display & Brightness settings.
- 9 On the right side, scroll down and tap **View**.



Note

This setting is not available on iPads smaller than 11 inches.

- 10 Select **Zoomed** or **Larger Text**, then tap **Set**.
- 11 In the popup, tap **Use Zoomed**. The screen will reset briefly.
- 12 Select **Cancel** to close the Use Zoomed popup.
- 13 You must now do a full iPad restart to apply the Zoom settings. Tap **General** on the left, then on the right swipe down and tap **Shut Down**. When the iPad has powered off, press the Power button to turn it on again.



What does AssistiveTouch have to do with eye gaze?

AssistiveTouch is designed for people who have difficulty touching the screen. The AssistiveTouch menu allows you to do “touch” functions, such as tap and scroll, using eye gaze. It also provides gaze-accessible shortcuts to things like the Home screen and App Switcher, which are typically accessed through gestures.

Set up AssistiveTouch



- 1 Tap **Settings**.



- 2 On the left side, tap **Accessibility**.
- 3 On the right side, tap **Touch**.
- 4 Tap **AssistiveTouch**, then turn it on.

Eye gaze is now enabled. You will see the Pointer, which shows the location of your eye gaze. The AssistiveTouch menu button also appears on the screen.




Note

If you are prompted to customize the Top Level Menu, tap **Cancel**.

- 5 Drag the AssistiveTouch menu button to the top third of the screen on the right side.

Customize the AssistiveTouch Menu

- 1 On the right side, select **Customize Top Level Menu**.
- 2 Tap **+** to change the number of icons to 8.
- 3 Tap **Notification Center**.
- 4 Swipe to the bottom of the list and tap **Toggle Pause/Resume Dwell**. Tap outside the list to close it.
- 5 Tap the **Gestures** button.
- 6 Swipe to and then tap **Move Menu**. Tap anywhere outside the menu to close it.
- 7 Continue editing the menu icons until your menu matches the one shown to the right.
- 8 On the left, select **Accessibility**.
- 9 On the right, select **Touch**.
- 10 On the right, select **AssistiveTouch**.



Configure Dwell Control

- 1 On the right side, swipe to the bottom of the AssistiveTouch menu. Turn on **Dwell Control**.



Note

If you are prompted to customize the AssistiveTouch Top Level Menu, tap **No**.

- 2 Tap the minus sign (–) next to **Seconds** to change the dwell time to 1.5 seconds.



Tip

This dwell time setting is to get you started. You can change the dwell time again later to suit your needs.

Waking and Unlocking Your iPad


The iPad wakes automatically when the eye tracker detects the user's eyes. Unlock the iPad by selecting the AssistiveTouch Button, then Home.

Try it now:

- 1 Press the top button to lock the iPad.
- 2 Look at the screen for a few moments. The iPad wakes and displays the lock screen.
- 3 Tap or use your eyes to select the **AssistiveTouch menu** button then select **Home**.



The iPad is now unlocked.



Tip

If you are interested in extra security on your device, we recommend using Face ID. Face ID allows you to unlock your device without entering a passcode each time. Set up Face ID in *iPad OS Settings > Face ID & Passcode*.

Set Up TD Snap


- 1 Select the **TD Snap** app.




- 2 Follow the prompts to set up your TD Snap user. Remember to choose the Eye Gaze access method.

Mount and Position


Prepare the user to start using eye tracking by positioning them comfortably. If they use glasses, make sure that they are wearing them and that the lenses are clean.



Place the TD Pilot on your mounting system or a stable surface in front of the user at, or slightly below, their eye level. If their head is tilted left or right, tilt the TD Pilot to match. It is important that the angle of the screen surface match the angle of the user's face. Most users, when seated at a table or desk, will need to have the TD Pilot positioned higher than the table surface.



You may need to refine the position of the device during the Calibrate steps on the next page. Always adjust the position of the device to suit the user, not the other way around.



Note

A mounting system is the best option for precise device positioning that is easy to adjust throughout the day. Several mounting options are available, including floor mounts, desk mounts, and wheelchair mounts. Visit [TobiiDynavox.com](https://www.tobii-dynavox.com) or contact your local Tobii Dynavox partner.

Calibrate

- 1 Tap the **CoPilot** app.



- 2 Follow the on-screen instructions in CoPilot to calibrate the eye tracker to the user's eyes.
- 3 When the user has finished calibrating, swipe up from the bottom of the screen to return to the Home screen.

Ready!



Your TD Pilot is now ready to use. Launch TD Snap to start communicating.

Continue to the TD Snap Training Cards

Continue your set up journey with the TD Snap Training Cards that came in the box with your TD Pilot.

The training cards walk you through the features of TD Snap, basic editing, backing up your data, and some ideas to help you integrate TD Snap into your daily life.



Tip

Customize your TD Snap eye gaze settings in *TD Snap > Edit > User > Access Method*.

Additional Resources

Scan the QR codes or use the links.



myTobiiDynavox
mytobiidynavox.com



Apple
apple.com/accessibility



TD Facebook Community
qrco.de/TDFB



Tobii Dynavox Learning Hub
learn.tobiidynavox.com



TD Pilot User's Manual
qrco.de/PilotDocs



TD Pilot Support Page
qrco.de/PilotHelpUK



Documentation
TD CoPilot > Settings > Help



UK Technical Support
0114 481 0011
support.uk@tobiidynavox.com