

Release Notes Boardmaker Studio

Release Version 1.5

Release date: 2018-10-22

New Features

- Sync activities between your local computer and Boardmakeronline.com — you can work offline on your local machine, then sync changes to Boardmakeronline.com the next time you have an internet connection.
- Upload activities to Boardmakeronline.com right from the Boardmaker Studio software.
- New Dashboard UI (formerly “Connect”) makes it easier to find, edit, and play activities.
- Full single sign on integration — use the same log in information that you use on other Tobii Dynavox and Boardmaker websites.

Improvements

- Rebranding.
- Refreshed UI.
- Bug fixes.

Release Version 1.4

Release date: 2015-01-22

New Features

- The new Symbol Editor allows users to modify symbol images to better customize their projects.
- Main user interface has been refreshed to provide consistency between Boardmaker Studio and Boardmaker Online Editor.
- Added option to choose Portrait or Landscape orientation from New Page dialog screen.
- New Message Window option allows content to use global Symbolate setting to determine if symbols should be used.

Improvements

- Bug fixes.

Release Version 1.3

Release date: 2014-05-08

New Features

- The Connect screen has been updated to integrate Boardmaker Studio with *BoardmakerAchieve.com*.
- A button has been added to the *Connect* screen (**Convert v5/v6 Boards**) to quickly convert a Boardmaker v5 or v6 board or set of boards (*.bm2*, *.zbp*, or *.zip*) into a Boardmaker Studio project file (*.bpf*) directly from the *Connect* screen.
- Screen objects in Boardmaker Studio can now be copied to the Windows clipboard.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- Several sets of templates (Question and Answer, Word Study, and Sorting) have been modified. These "performance enabled templates" can track student performance and report back results when the activities are used by students in a Boardmaker Achieve subscription account. These templates are indicated by a red "P" icon.
- When creating a new activity from a template, it is no longer necessary to name the activity prior to creating it.

Improvements

- Several actions and functions have been added to the *Screen Objects* category in the Action Editor.
- Bug fixes.

Release Version 1.2

Release date: 2012-01-10

New Features

- This version of the software includes a new gadget category, *Quick Edit* gadgets. Gadgets in this category allow the user to quickly modify the text and symbols on buttons while in Play mode by copying text or symbols from one button to another "on-the-fly."
- Two new categories of templates, *Classroom Routines* and *Graphic Organizers*, have been added to this version of the software. New templates have also been added to the *Word Study* and *Sorting* categories and to the *Mathematics/Money* folder in the U.K., Canadian, Australian, and New Zealand versions.

Improvements

- New gadgets have been added to the *Draggable Buttons/Sorting* category and the *Fun Interactions* category.
- New templates have been added to the *Word Study* folder (U.S. and Canada only), *Sorting*, and *Mathematics/Money* (U.K., Canada, Australia, and New Zealand only).
- Bug fixes.

Release Version 1.1

Release date: 2010-12-17

New Features

- New Access Method: 2-switch scanning.
- New object property: Target Scanning. This option controls the size of the positioning grid that appears when placing draggable objects inside a target using the Scanning access method.
- Added audio cue voice setting, so that audio cues can be spoken in a different voice than speaking voice.

Improvements

- Draggable objects are now accessible in Play mode through any Access Method.
- A *Release Time* slider bar has been added to the Touch Exit settings.
- New Gadgets and Templates.
- New Actions and Functions.
- Bug fixes and enhancements.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact