

Release Notes Tobii Dynavox Hardware Listener

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 1.8.2

Release date: 2024-02-27

Improvements

- Fixed an issue for I-Series devices where the Adaptive Button action “Enable/Disable Eye Tracker Pause” would not work in TD Snap. Requires TD Snap version 1.31 or above.

Release Version 1.8.0

Release date: 2023-03-06

Improvements

- Fixed an issue where logging was not working.

Release Version 1.7.0

Release date: 2022-06-13

New Features

- Verifies Partner Window driver at startup and installs correct one if not present.

Release Version 1.6.0

Release date: 2022-02-01

Improvements

- Resolved issue with TD Snap not restoring after closing Track Status.

Release Version 1.5.4

Release date: 2021-12-07

Improvements

- Bug fixes

Release Version 1.5.3

Release date: 2021-07-06

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: <https://www.tobiidynavox.com/pages/product-support>

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit <https://www.tobiidynavox.com/pages/contact-us>.

Improvements

- Bug fixes

Release Version 1.5.2

Release date: 2021-05-10

New Features

- Automatic headphone jack detection (I-13/I-16 only)
If enabled, the audio output device is automatically changed when headphones are plugged in or unplugged.

Release Version 1.4.0

Release date: 2021-02-08

New Features

- Ability to assign mouse click actions to switches (I-13/I-16 only)
 - Left Click
 - Middle Click
 - Right Click
- Ability to assign start program to switches (I-13/I-16 only)
 - Select a program to start when a switch is selected
- Ability to assign volume actions to switches (I-13/I-16 only)
 - Enabled/Disable Mute
 - Volume Down
 - Volume Up
- Ability to assign eye tracker actions to switches (I-13/I-16 only)
 - Show/Hide Track Status
 - Enable/Disable Eye Tracker Pause
 - Calibrate Eye Tracker
 - Enable/Disable Computer Control

Release Version 1.3.4

Release date: 2020-12-07

Improvements

- Bug fixes

Release Version 1.3.3

Release date: 2020-11-10

New Features

- Ability to assign power off device action to switches (I-13/I-16 only)
 - Press and hold 2-7 seconds for normal Windows shutdown

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product.
Find our Support pages online at: <https://www.tobiidynavox.com/pages/product-support>

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit <https://www.tobiidynavox.com/pages/contact-us>.

- Press and hold 7 seconds or more for forced Windows shutdown

Improvements

- Bug fixes

Release Version 1.3.1

Release date: 2020-08-10

Improvements

- Bug fixes

Release Version 1.3.0

Release date: 2020-06-08

New Features

- Ability to assign keyboard key actions to switches (I-13/I-16 only)

Release Version 1.2.0

Release date: 2020-04-08

Improvements

- Hardware Listener separated from Hardware Settings into its own installer

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product.
Find our Support pages online at: <https://www.tobiidynavox.com/pages/product-support>

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit <https://www.tobiidynavox.com/pages/contact-us>.