

Release Notes Tobii Dynavox Hardware Listener

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 1.3.4.243

Release date: 2020-12-07

Improvements

- Bug fixes

Release Version 1.3.3.240

Release date: 2020-11-10

New Features

- Ability to assign power off device action to switches (I-13/I-16 only)
 - Press and hold 2-7 seconds for normal Windows shutdown
 - Press and hold 7 seconds or more for forced Windows shutdown

Improvements

- Bug fixes

Release Version 1.3.1.234

Release date: 2020-08-10

Improvements

- Bug fixes

Release Version 1.3.0.230

Release date: 2020-06-08

New Features

- Ability to assign keyboard key actions to switches (I-13/I-16 only)

Release Version 1.2.0.223

Release date: 2020-04-08

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Improvements

- Hardware Listener separated from Hardware Settings into its own installer

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact