

## Release Notes Tobii Dynavox Hardware Settings

**Customer support:** For more information or support, please visit [www.TobiiDynavox.com](http://www.TobiiDynavox.com).

### Release Version 1.3.4.243

**Release date:** 2020-12-07

#### Improvements

- Updates to Snap log collection
- Add PODD and Gateway to Product Key information

### Release Version 1.3.3.240

**Release date:** 2020-11-10

#### New Features

- Ability to assign power off device action to switches (I-13/I-16 only)
  - Press and hold 2-7 seconds for normal Windows shutdown
  - Press and hold 7 seconds or more for forced Windows shutdown

#### Improvements

- Keyboard key translations fixes
- Improved log collection
- Bug fixes

### Release Version 1.3.1.234

**Release date:** 2020-08-10

#### Improvements

- Add keyboard key translations
- Collect Configuration Manager logs

### Release Version 1.3.0.230

**Release date:** 2020-06-08

#### New Features

- Ability to assign keyboard key actions to switches (I-13/I-16 only)

#### Support for Your Tobii Dynavox Device

##### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: [www.TobiiDynavox.com/support-training](http://www.TobiiDynavox.com/support-training)

##### Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)

## Improvements

- Improved log collection
- Bug fixes

## Release Version 1.2.0.223

**Release date:** 2020-04-08

## New Features

- Ability to run on EM12, I-110, Indi, and Indi7

## Improvements

- Display Eye Tracker firmware version in System & Diagnostics
- Update default Square Button setting to “Toggle Pause”
- Hardware Listener removed and separated into its own installer

## Release Version 1.1.0.205

**Release date:** 2019-12-16

## New Features

- Ability to disable the use of touch on the adaptive volume buttons.

## Improvements

- Improved log collection
- Bug fixes

## Release Version 1.0.1.196

**Release date:** 2019-10-18

## Improvements

- UI string updates to Swedish
- Bug fixes

## Support for Your Tobii Dynavox Device

### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: [www.TobiiDynavox.com/support-training](http://www.TobiiDynavox.com/support-training)

### Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)