

Release Notes Tobii Dynavox Hardware Settings

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 1.4.0.249

Release date: 2021-02-08

New Features

- Ability to assign mouse click actions to switches (I-13/I-16 only)
 - Left Click
 - Middle Click
 - Right Click
- Ability to assign start program to switches (I-13/I-16 only)
 - Select a program to start when a switch is selected
- Ability to assign volume actions to switches (I-13/I-16 only)
 - Enabled/Disable Mute
 - Volume Down
 - Volume Up
- Ability to assign eye tracker actions to switches (I-13/I-16 only)
 - Show/Hide Track Status
 - Enable/Disable Eye Tracker Pause
 - Calibrate Eye Tracker
 - Enable/Disable Computer Control

Release Version 1.3.4.243

Release date: 2020-12-07

Improvements

- Updates to Snap log collection
- Add PODD and Gateway to Product Key information

Release Version 1.3.3.240

Release date: 2020-11-10

New Features

- Ability to assign power off device action to switches (I-13/I-16 only)
 - Press and hold 2-7 seconds for normal Windows shutdown
 - Press and hold 7 seconds or more for forced Windows shutdown

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Improvements

- Keyboard key translations fixes
- Improved log collection
- Bug fixes

Release Version 1.3.1.234

Release date: 2020-08-10

Improvements

- Add keyboard key translations
- Collect Configuration Manager logs

Release Version 1.3.0.230

Release date: 2020-06-08

New Features

- Ability to assign keyboard key actions to switches (I-13/I-16 only)

Improvements

- Improved log collection
- Bug fixes

Release Version 1.2.0.223

Release date: 2020-04-08

New Features

- Ability to run on EM12, I-110, Indi, and Indi7

Improvements

- Display Eye Tracker firmware version in System & Diagnostics
- Update default Square Button setting to "Toggle Pause"
- Hardware Listener removed and separated into its own installer

Release Version 1.1.0.205

Release date: 2019-12-16

New Features

- Ability to disable the use of touch on the adaptive volume buttons.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Improvements

- Improved log collection
- Bug fixes

Release Version 1.0.1.196

Release date: 2019-10-18

Improvements

- UI string updates to Swedish
- Bug fixes

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact