

## Release Notes Tobii Dynavox PCEye 5 Installer

### Introduction

The Tobii Dynavox PCEye 5 Installer is an installation bundle for the software needed for running the Tobii Dynavox PCEye 5 eye tracker.

### System Requirements

- Windows 10
- .NET Framework 4.5 Full

**Customer support:** For more information or support, please visit [www.TobiiDynavox.com](http://www.TobiiDynavox.com).

### Release Version 1.0.11

**Release date:** 2023-08-15

#### Improvements

- Changed bundle installer name.
- Updated Switcher version to 2.2.1.2423.
- Updated Tobii Experience Software for Windows (PCEye5) version to 4.180.0.29190.
- Updated TD Control version to 1.79.1.837.
- Updated Tobii Dynavox Eye Tracking version to 3.1.7.2703.

### Release Version 1.0.10

**Release date:** 2023-02-21

#### Improvements

- Updated error message that is shown when installation failed.
- Updated Switcher version to 2.1.0.62303.
- Updated Update Notifier version to 1.8.1.62923.
- Updated TD Control version to 1.75.0.62106.
- Updated Tobii Dynavox Eye Tracking version to 3.1.4.62985.

### Release Version 1.0.9

**Release date:** 2022-11-03

#### Improvements

- Updated Display Setup configuration step to notify the user if eye tracker is not connected.
- Updated Switcher version to 2.0.0.58951.
- Updated Tobii Experience Software for Windows (PCEye5) version to 4.163.0.24281.
- Updated Tobii Dynavox Eye Tracking to 3.1.3.58139.
- Updated TD Control version to 1.69.0.58518.

### Release Version 1.0.8

**Release date:** 2022-05-24

## Improvements

- Updated Switcher version to 1.9.0.53636.
- Updated Tobii Experience Software for Windows (PCEye5) version to 4.140.0.19139.
- Updated Update Notifier version to 1.8.0.52281.
- Updated Tobii Dynavox Eye Tracking to 3.1.2.54134.
- Updated TD Control version to 1.59.0.53936.

## Release Version 1.0.7

**Release date:** 2022-01-31

## Improvements

- Fixed inconsistency in the Swedish translations.
- Updated Switcher version to 1.7.0.47736.
- Updated Tobii Experience Software for Windows (PCEye5) version to 4.127.0.16441.
- Updated Update Notifier version to 1.7.4.501.
- Updated TD Control version to 1.51.0.48512.

## Release Version 1.0.6

**Release date:** 2021-10-05

## Improvements

- Updated some text for Chinese and Hebrew languages.
- Updated Computer Control in the installer to 1.43.0.42495 version.
- Updated Tobii Dynavox Eye Tracking in the installer to 3.1.0.40640 version.
- Updated Tobii Dynavox Update Notifier in the installer to 1.7.3.488 version.
- Added Switcher 1.2.0.43071 to the installer.

## Release Version 1.0.5

**Release date:** 2021-06-29

## Improvements

- Updated Computer Control in the installer to 1.37.0.37801 version.
- Updated Tobii Dynavox Eye Tracking in the installer to 3.0.9.367 version.
- Updated Tobii Dynavox Update Notifier in the installer to 1.7.2.481 version.
- Updated Tobii Experience Software for Windows (PCEye5) in the installer to 4.107.0.12798 version.

## Release Version 1.0.3

**Release date:** 2020-11-30

## Improvements

- Updated Computer Control in the installer to 1.23.0.27322 version.
- Updated Tobii Dynavox Eye Tracking in the installer to 3.0.8.313 version.

- Updated Tobii Experience Software for Windows (PCEye5) in the installer to 4.70.0.7919 version

## Support for Your Tobii Dynavox Device

### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product.  
Find our Support pages online at: [www.TobiiDynavox.com/support-training](http://www.TobiiDynavox.com/support-training)

### Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)