


The Boardmaker Player in Snap – Adding Activities: Quick Guide

The integrated Boardmaker Player in Snap allows you to play preloaded and interactive activities from Boardmaker right in your Snap software. Preloaded activities include an accessible bookshelf, accessible games, an interactive calendar, math tools, as well as samples from Core First Learning, Reading Avenue, and Expedition Education curriculum titles. This document will explain how to add activities from your Boardmaker account to be used in the integrated Boardmaker Player in Snap. Your accounts are connected by your single sign on email address; therefore, you must be logged into myTobiiDynavox in Snap using the same email address to log into your Boardmaker account.


Note: Activities played in the integrated Boardmaker Player in Snap do not connect to the Boardmaker Student Center and do not collect data.



e. Log in using your Tobii Dynavox single sign on email address and password

Note: Use the same email address that is logged into myTobiiDynavox in Snap in order to connect the two accounts

f. Select Log In



2. Add activities from the Community to My Boardmaker


Search the Community

- a. Click on the **magnifying glass** in the menu bar
- b. Type in your search term (i.e.: activity type, topic, etc.)
- c. Hit enter on your keyboard
- d. Scroll through or user filters to find activities
- e. Select **the Title** of the activity to see a preview

 Select the + to add the activity to My Boardmaker

Note: The blue tool bar will change when the activity is successfully added to My Boardmaker

g. Once successfully added, select the **My Boardmaker** tab



3. Add activities from My Boardmaker to the integrated Boardmaker Player


Add any activity to the Snap Activities Folder

- a. From the My Boardmaker tab, locate the activity you want to add
- b. Select the **Move activity to folder** icon

c. From the popup select **Snap Activities** or an appropriate subfolder within Snap Activities

- 4. Open Snap
- 5. Navigate to the Boardmaker Library
- If connected to Wi-Fi and logged into the associated myTobiiDynavox account in Snap your activity is now available to be used in the integrated Boardmaker Player

Note: If the Boardmaker Library was already open in Snap you may need to navigate out and reopen to see the addition of the activity



Note: Any activity added to the Snap Activity folder in My Boardmaker will appear in the integrated Boardmaker Player of any device that is logged into the mytobiidynavox account using the same email address. If you are working with multiple AAC users, want to assign specific activities to each user or want to collect data click https://get.boardmakeronline.com to learn how to get a free trial or subscribe to Boardmaker.