

Release Notes Tobii Dynavox Eye Tracking Settings

Introduction

Tobii Dynavox Eye Tracking is a software solution for setting up and configuring Eye tracking on the I-Series I-13/I-16 devices and also for the PCEye 5. The release notes contain information for I-Series I-13/I-16 and for PCEye 5.

System Requirements

- Windows 10
- .NET Framework 4 Full
- Hardware
 - I-Series I-13/I-16
 - PCEye 5

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 3.1.3

Release date: 2022-10-11

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Made sure that the pause indicator is not shown on top of certain Tobii Dynavox applications.
- Added a possibility for client applications to get information about eye tracker connection state.
- Fixed a possible crash when using eye gaze to click on buttons.

Release Version 3.1.2

Release date: 2022-05-24

Release channels:

- Tobii Dynavox Update Notifier

New Features

- Added a possibility to enable or disable sound feedback during calibration process.

Improvements

- Made sure that eye tracking pause indicator is not shown when certain Tobii Dynavox applications are running.
- Fixed Arabic application text.
- Fixed Japanese application text.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Release Version 3.1.1

Release date: 2022-02-15

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Updated display setup implementation and added support for PCEye Brackets
- Made sure that TD Control is hibernated whenever TD Eye Tracking Settings application is shown
- Added sound feedback to the eye tracking calibration process
- Updated the list of log files that are collected by the application for further troubleshooting
- Fixed a bug where enabling continuous dwell activation setting externally didn't update the application user interface
- Updated translations for Japanese language

Release Version 3.1.0

Release date: 2021-08-17

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Eye tracking off-screen menu is shown depending on other Tobii Dynavox applications running at the same time.

Release Version 3.0.9

Release date: 2021-03-03

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Made the eye tracking pause indicator more visible
- Disabled showing eye tracking off screen menu when other Tobii Dynavox applications are running
- Fixed a bug where the application processes were not started automatically after the upgrade

Release Version 3.0.8

Release date: 2020-11-30

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Added a possibility to change opacity of the region interaction feedback

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

- Added a possibility to get notifications whenever missing eye tracking components are not installed
- Added a possibility to generate troubleshooting logs using eye gaze
- Fixed a bug where sometimes creating default calibration set eye tracker in the power saving mode
- Fixed a bug where sometimes eye gaze did not work when connecting eye tracker after startup of the application

Release Version 3.0.7

Release date: 2020-09-23

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Added support for metrics
- Added a prevention from installing the software if Tobii Gaze Interaction Software, Gaze Point 2 or Windows Control 2 is installed on the device
- Removed Wake-On-Gaze option in settings for peripheral users
- Fixed a bug where default calibration profile was created before display setup has been done
- Fixed a bug where display setup option in settings was enabled when the tracker was not connected
- Disabled pop up messages from the software
- Added a possibility to generate troubleshooting files from the software
- Added translation to Arabic and Russian languages

Release Version 3.0.6

Release date: 2020-03-30

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Added a possibility to open Computer Control Review version from off-screen menu
- Added a possibility to set up eye tracking display configuration
- Updated switch implementation to activate on press event instead of release event
- Fixed a bug where the application was not properly updated through Update Notifier application
- Fixed incorrect Italian translations in Settings window
- Fixed a bug where eye tracker was disabled sometimes after a new profile with default calibration was created

Release Version 3.0.5

Release date: 2020-01-27

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Added a possibility to create a user calibration profile with default calibration

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact

Release Version 3.0.4

Release date: 2019-12-17

Release channels:

- Tobii Dynavox Update Notifier

Improvements

- Pop-up menu for Wake-On-Gaze/Sleep-On-Gaze to stay open longer
- Removed the text follow the targets it's OK to blink when starting a new calibration

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com/support-training

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Dynavox solution consultant or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact